



Product description

| | |
|---------------------|---|
| Order number | 62.732.519 |
| Product description | Screw cap tube, working volume: 15 ml, (LxØ): 76 x 20 mm, material: PP, skirted conical base, transparent, screw cap, natural, cap enclosed, label/print: white, with graduations, 100 piece(s)/bag |

Product characteristics

| | |
|--------------------------|----------------------|
| Colour of print/label | white |
| Graduation | yes |
| Closure type | screw cap |
| Base shape | skirted conical base |
| Volume of work | 15 ml |
| Centrifugation max (RCF) | 4500 x g |

Size

| | |
|-------------------|-------|
| Diameter | 20 mm |
| Length of product | 76 mm |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Material & colours

| | |
|-------------------|-----------------------------------|
| Product material | Polypropylene (PP) |
| Colour of product | transparent |
| Closure material | High Density Polyethylene (HD-PE) |
| Colour of cap | natural |

Purity & certification

| | |
|---------------------------|-------------------------|
| Satisfies the requirement | ADR |
| Product category | no medical device IVD |
| Batched | yes |

Packaging

| | |
|-------------------------------|---------------|
| Minimum order qty. | 500 |
| Type of smallest subpackaging | bag |
| Piece(s) / inner box | 100 |
| Piece (s) / outer case | 500 |
| Piece(s) / pallet | 18000 |
| Depth of case | 398 mm |
| Width of case | 398 mm |
| Height of case | 265 mm |
| Case volume | 0.042 cbm |
| Weight of product | 0.0064 kg |
| Weight of case | 3.81 kg |
| EAN of inner box | 4038917251617 |
| EAN case | 4038917251600 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

