



### Product description

Order number	62.9924.293
Product description	Tube, working volume: 10 ml, (LxØ): 100 x 16 mm, material: PP, conical base, transparent, screw cap, natural, cap enclosed, with print, label/print: white, with graduations, 500 piece(s)/bag

### Product characteristics

Label/ Print	with print
Colour of print/label	white
Graduation	yes
Closure type	screw cap
Base shape	conical base
Volume of work	10 ml
Centrifugation max (RCF)	4500 x g
Autoclavable (max. °C)	121 °C

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





### Size

Diameter	16 mm
Length of product	100 mm

### Material & colours

Product material	Polypropylene (PP)
Colour of product	transparent
Closure material	High Density Polyethylene (HD-PE)
Colour of cap	natural

### Purity & certification

Satisfies the requirement	IATA
Product category	no medical device   IVD
Batched	yes

### Packaging

Minimum order qty.	1000
Type of smallest subpackaging	bag
Piece(s) / inner box	500
Piece (s) / outer case	1000
Piece(s) / pallet	32000
Case volume	0.0559 cbm
Weight of product	0.0005 kg
EAN of inner box	4038917369923
EAN case	4038917534321

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

