



**Product Specification**  
**Screw cap tube, 25 ml, (LxØ): 90 x 25 mm, PP,**  
**with paper label**



**Product description**

Order number	63.9922.252
Product description	Screw cap tube, working volume: 25 ml, (LxØ): 90 x 25 mm, material: PP, skirted conical base, transparent, screw cap, white, cap assembled, with paper label, label/print: white/black, Patient information, with graduations, 500 piece(s)/bag

**Product characteristics**

Label/ Print	with paper label
Colour of print/label	white/black
Graduation	yes
Closure type	screw cap
Base shape	skirted conical base
Volume of work	25 ml
Centrifugation max (RCF)	4500 x g

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Screw cap tube, 25 ml, (LxØ): 90 x 25 mm, PP,  
with paper label**

**Page 2**

**Size**

Diameter 25 mm  
Length of product 90 mm

**Material & colours**

Product material Polypropylene (PP)  
Colour of product transparent  
Closure material High Density Polyethylene (HD-PE)  
Colour of cap white

**Purity & certification**

Satisfies the requirement IATA  
Product category In vitrodiagnostic device, CE  
CE certified CE  
Batched yes

**Packaging**

Minimum order qty. 500  
Type of smallest subpackaging bag  
Piece(s) / inner box 500  
Piece (s) / outer case 500  
Piece(s) / pallet 15000  
Depth of case 608 mm  
Width of case 324 mm  
Height of case 264 mm  
Case volume 0.052 cbm  
Weight of product 0.0084 kg  
Weight of case 4.98 kg  
EAN of inner box 4038917151559  
EAN case 4038917027779

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-07-03

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**