



Product Specification

Multi-purpose container, max. volume: 250 ml,
(LxØ): 78 x 70 mm, graduated, PP, transparent,
with safety label

Page 1



Product description

Order number	75.9922.518
Product description	Multi-purpose container, max. volume: 250 ml, (LxØ): 78 x 70 mm, transparent, cap: yellow, graduated, material: PP, screw cap, cap assembled, with safety label, sterile, 240 piece(s)/bag

Product characteristics

Label/ Print	with safety label
Graduation	yes
Closure type	screw cap
Base shape	flat base
Volume of work	250 ml

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-08-18

This document was prepared by EDP support and is valid without signature.



SARSTEDT



Product Specification

**Multi-purpose container, max. volume: 250 ml,
(LxØ): 78 x 70 mm, graduated, PP, transparent,
with safety label**

Page 2

Size

Sample volume	250 ml
Max. volume	250 ml
Diameter	70 mm
Height of product	78 mm
Length of product	78 mm

Material & colours

Product material	Polypropylene (PP)
Colour of product	transparent
Closure material	High Density Polyethylene (HD-PE)
Colour of cap	yellow

Purity & certification

Product category	no medical device IVD
Purity standard	sterile
Sterilisation	Gamma irradiation
Batched	yes

Packaging

Minimum order qty.	240
Type of smallest subpackaging	bag
Piece(s) / inner box	240
Piece (s) / outer case	240
Piece(s) / pallet	1440
Depth of case	646 mm
Width of case	450 mm
Height of case	498 mm
Case volume	0.1448 cbm
Weight of product	0.0404 kg
Weight of case	10.93 kg
EAN of inner box	4038917250955
EAN case	4038917350914

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-08-18

This document was prepared by EDP support and is valid without signature.



SARSTEDT