



### **Product description**

Order number	86.1170.300
Product description	Transfer pipette, nominal volume: 1 ml, (LxW): 115 x 10 mm, material: LD-PE, transparent, droplet size: 40 µl - 50 µl, 500 piece(s)/case

### **Product characteristics**

Graduation	yes
------------	-----

### **Size**

Sample volume	1 ml
Width of product	10 mm
Length of product	115 mm
Drop size	40 µl - 50 µl

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





### **Material & colours**

Product material Low Density Polyethylene (LD-PE)  
Colour of product transparent

### **Purity & certification**

Product category no medical device | IVD  
Batched yes

### **Packaging**

Minimum order qty. 500  
Type of smallest case  
subpackaging  
Piece(s) / inner box 500  
Piece (s) / outer case 4000  
Piece(s) / pallet 156000  
Case volume 0.0401 cbm  
Weight of product 0.001 kg  
Weight of case 4.55 kg  
EAN of inner box 4038917378352  
EAN case 4038917378345

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

