



Product description

Order number	86.1175.001
Product description	Transfer pipette, nominal volume: 6 ml, (LxW): 146 x 15 mm, material: LD-PE, transparent, droplet size: 20 μ l - 30 μ l, sterile, 1 piece(s)/blister

Product characteristics

Graduation no

Size

Sample volume	6 ml
Width of product	15 mm
Length of product	146 mm
Drop size	20 μ l - 30 μ l

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Material & colours

Product material Low Density Polyethylene (LD-PE)
Colour of product transparent

Purity & certification

Product category In vitrodiagnostic device, CE
CE certified CE
Purity standard sterile
Sterilisation Gas sterilisation
Batched yes

Packaging

Minimum order qty. 840
Type of smallest blister
subpackaging
Piece(s) / inner box 42
Piece (s) / outer case 840
Piece(s) / pallet 23520
Depth of case 533 mm
Width of case 333 mm
Height of case 250 mm
Case volume 0.0444 cbm
Weight of product 0.0016 kg
Weight of case 2.56 kg
EAN of inner box 4038917126311
EAN case 4038917020589

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

