




**Product Specification**  
**Serological pipette, with tip, plugged, 10 ml,  
sterile, 25 piece(s)/bag**

**Page 1**



**Product description**

Order number 86.1254.025  
Product description Serological pipette, with tip, plugged, 10 ml, coding colour: orange, material: PS, sterile, 25 piece(s)/bag

**Product characteristics**

Label/ Print with printed graduations  
Colour of print/label orange/black  
Graduation yes  
Design with tip, plugged  
Volume of work 1 ml - 13 ml  
Coding colour orange

**Size**

Sample volume 10 ml  
Diameter 10 mm  
Length of product 335 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Serological pipette, with tip, plugged, 10 ml,  
sterile, 25 piece(s)/bag**

**Page 2**

**Material & colours**

Product material Polystyrene (PS)  
Colour of product transparent

**Purity & certification**

Product category no medical device | IVD  
Purity standard sterile  
Sterilisation Gamma irradiation  
Batched yes

**Packaging**

Minimum order qty. 500  
Type of smallest bag  
subpackaging  
Piece(s) / inner box 25  
Piece (s) / outer case 500  
Piece(s) / pallet 27000  
Depth of case 413 mm  
Width of case 250 mm  
Height of case 270 mm  
Case volume 0.0279 cbm  
Weight of product 0.0091 kg  
Weight of case 5.26 kg  
EAN of inner box 4038917127318  
EAN case 4038917021135

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-07-17

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**