




**Product Specification**  
**Serological pipette, with tip, plugged, 25 ml,  
sterile, non-pyrogenic/endotoxin-free,  
non-cytotoxic, 1 piece(s)/blister**

**Page 1**



**Product description**

Order number	86.1685.001
Product description	Serological pipette, with tip, plugged, 25 ml, coding colour: red, material: PS, sterile, non-pyrogenic/endotoxin-free, non-cytotoxic, 1 piece(s)/blister

**Product characteristics**

Label/ Print	with printed graduations
Colour of print/label	red/black
Graduation	yes
Design	with tip, plugged
Volume of work	3 ml - 35 ml
Coding colour	red

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Serological pipette, with tip, plugged, 25 ml,  
sterile, non-pyrogenic/endotoxin-free,  
non-cytotoxic, 1 piece(s)/blister**

**Page 2**

**Size**

Sample volume	25 ml
Diameter	15 mm
Length of product	345 mm

**Material & colours**

Product material	Polystyrene (PS)
Colour of product	transparent

**Purity & certification**

Product category	no medical device   IVD
Purity standard	sterile, non-pyrogenic/endotoxin-free, non-cytotoxic
Sterilisation	Gamma irradiation
Pyrogen-free	yes
Batched	yes

**Packaging**

Minimum order qty.	200
Type of smallest subpackaging	blister
Piece(s) / inner box	25
Piece (s) / outer case	200
Piece(s) / pallet	10800
Depth of case	420 mm
Width of case	250 mm
Height of case	290 mm
Case volume	0.0303 cbm
Weight of product	0.0169 kg
Weight of case	3.98 kg
EAN of inner box	4038917133012
EAN case	4038917024655

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-07-18

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**