




**Product Specification**  
**Serological pipette, with tip, plugged, 25 ml,  
sterile, 20 piece(s)/bag**

**Page 1**



**Product description**

Order number 86.1685.020  
Product description Serological pipette, with tip, plugged, 25 ml, coding colour: red, material: PS, sterile, 20 piece(s)/bag

**Product characteristics**

Label/ Print with printed graduations  
Colour of print/label red/black  
Graduation yes  
Design with tip, plugged  
Volume of work 3 ml - 35 ml  
Coding colour red

**Size**

Sample volume 25 ml  
Diameter 15 mm  
Length of product 345 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Serological pipette, with tip, plugged, 25 ml,  
sterile, 20 piece(s)/bag**

**Page 2**

**Material & colours**

Product material Polystyrene (PS)  
Colour of product transparent

**Purity & certification**

Product category no medical device | IVD  
Purity standard sterile  
Sterilisation Gamma irradiation  
Batched yes

**Packaging**

Minimum order qty. 200  
Type of smallest bag  
subpackaging  
Piece(s) / inner box 20  
Piece (s) / outer case 200  
Piece(s) / pallet 10800  
Depth of case 413 mm  
Width of case 250 mm  
Height of case 270 mm  
Case volume 0.0279 cbm  
Weight of product 0.0152 kg  
Weight of case 3.61 kg  
EAN of inner box 4038917133029  
EAN case 4038917040785

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-07-17

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**