



Product description

Order number 65.716.323
Product description Screw cap, red, sterile, suitable for screw cap micro tubes, 500 piece(s)/double bag

Product characteristics

Cap flat cap
Closure type screw cap
Suitable for screw cap micro tubes

Material & colours

Product material Polypropylene (PP)
Colour of product red
O-ring material Ethylene propylene diene rubber (EPDM)

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Purity & certification

Product category	In vitro diagnostic device, CE
CE certified	CE
Purity standard	sterile
Sterilisation	Electron irradiation
Batched	yes

Packaging

Type of smallest subpackaging	double bag
Piece(s) / inner box	500
Piece (s) / outer case	5000
Piece(s) / pallet	540000
Depth of case	351 mm
Width of case	280 mm
Height of case	138 mm
Case volume	0.0136 cbm
Weight of case	2.05 kg
EAN of inner box	4038917923415
EAN case	4038917923422

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

