



Product description

| | |
|---------------------|---|
| Order number | 72.689.006 |
| Product description | Reaction tube, working volume: 2 ml, material: PP, transparent, cap: violet, retaining loop cap, cap attached, 500 piece(s)/bag |

Product characteristics

| | |
|--------------------------|--------------------|
| Graduation | no |
| Cap | retaining loop cap |
| Volume of work | 2 ml |
| Centrifugation max (RCF) | 20000 x g |

Size

| | |
|-------------------|---------|
| Diameter | 10.8 mm |
| Length of product | 39 mm |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Material & colours

| | |
|-------------------|----------------------------------|
| Product material | Polypropylene (PP) |
| Colour of product | transparent |
| Closure material | Low Density Polyethylene (LD-PE) |
| Colour of cap | violet |

Purity & certification

| | |
|------------------|-------------------------------|
| Product category | In vitrodiagnostic device, CE |
| CE certified | CE |
| Batched | yes |

Packaging

| | |
|-------------------------------|---------------|
| Minimum order qty. | 5000 |
| Type of smallest subpackaging | bag |
| Piece(s) / inner box | 500 |
| Piece (s) / outer case | 5000 |
| Piece(s) / pallet | 150000 |
| Depth of case | 398 mm |
| Width of case | 398 mm |
| Height of case | 340 mm |
| Case volume | 0.0539 cbm |
| Weight of product | 0.0014 kg |
| Weight of case | 7.48 kg |
| EAN of inner box | 4038917116466 |
| EAN case | 4038917015158 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

