



### **Product description**

Order number	72.690.300
Product description	Reaction tube, working volume: 1.5 ml, material: PP, transparent, cap: transparent, retaining loop cap, cap attached, with injected graduation and writing space, 500 piece(s)/bag

### **Product characteristics**

Graduation	yes
Cap	retaining loop cap
Volume of work	1,5 ml
Centrifugation max (RCF)	20000 x g
Autoclavable (max. °C)	121 °C

### **Size**

Sample volume	1.5 ml
Length of product	39 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





### **Material & colours**

Product material	Polypropylene (PP)
Colour of product	transparent
Colour of cap	transparent

### **Purity & certification**

Product category	no medical device   IVD
Batched	yes

### **Packaging**

Minimum order qty.	500
Type of smallest subpackaging	bag
Piece(s) / inner box	500
Piece (s) / outer case	5000
Piece(s) / pallet	160000
Case volume	0.0558 cbm
Weight of product	0.0009 kg
Weight of case	5.14 kg
EAN of inner box	4038917378741
EAN case	4038917378734

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

