



### Product description

Order number	72.690.999
Product description	Reaction tube, working volume: 1.5 ml, material: PP, colour mix, cap: colour mix, retaining loop cap, cap attached, with injected graduation and writing space, 500 piece(s)/bag

### Product characteristics

Graduation	yes
Cap	retaining loop cap
Volume of work	1,5 ml
Centrifugation max (RCF)	20000 x g
Autoclavable (max. °C)	121 °C

### Size

Sample volume	1.5 ml
Length of product	39 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Material & colours**

Product material	Polypropylene (PP)
Colour of product	colour mix
Colour of cap	colour mix

**Packaging**

Minimum order qty.	1000
Type of smallest subpackaging	bag
Piece(s) / inner box	100
Piece (s) / outer case	5000
Piece(s) / pallet	140000
Depth of case	410 mm
Width of case	461 mm
Height of case	345 mm
Case volume	0.0649 cbm
Weight of case	5.19 kg
EAN of inner box	4038917116756
EAN case	4038917015288

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

