



### **Product description**

Order number	72.692.210
Product description	Screw cap micro tube, working volume: 1.5 ml, conical base, with knurling, transparent, cap: natural, cap attached, Biosphere® plus, 25 piece(s)/bag

### **Product characteristics**

Graduation	no
Cap	retaining loop cap
Closure type	screw cap
Knurls	yes
Base shape	conical base
Volume of work	1,5 ml
Centrifugation max (RCF)	20000 x g

### **Size**

Diameter	10.8 mm
Length of product	44 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





### **Material & colours**

Product material	Polypropylene (PP)
Colour of product	transparent
Closure material	Polypropylene (PP)
Colour of cap	natural
O-ring material	Ethylene propylene diene rubber (EPDM)

### **Purity & certification**

Satisfies the requirement	IATA
Product category	In vitro diagnostic device, CE
CE certified	CE
Purity standard	Biosphere® plus
Sterilisation	Gas sterilisation
Pyrogen-free	yes
Batched	yes

### **Packaging**

Minimum order qty.	500
Type of smallest subpackaging	bag
Piece(s) / inner box	250
Piece (s) / outer case	500
Piece(s) / pallet	55000
Depth of box	211 mm
Width of box	191 mm
Height of box	141 mm
Depth of case	398 mm
Width of case	223 mm
Height of case	160 mm
Case volume	0.0142 cbm
Weight of product	0.0015 kg
Weight of case	1.34 kg
EAN of inner box	4038917223775
EAN case	4038917223768

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

