



### **Product description**

Order number	72.733.001
Product description	Screw cap micro tube, working volume: 0.5 ml, conical base, with knurling, transparent, cap: natural, cap enclosed, 500 piece(s)/bag

### **Product characteristics**

Graduation	no
Cap	standard screw cap
Knurls	yes
Base shape	conical base
Volume of work	0,5 ml
Centrifugation max (RCF)	20000 x g

### **Size**

Diameter	8 mm
Length of product	44 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





### **Material & colours**

Product material	Polypropylene (PP)
Colour of product	transparent
Closure material	Polypropylene (PP)
Colour of cap	natural
O-ring material	Ethylene propylene diene rubber (EPDM)

### **Purity & certification**

Satisfies the requirement	IATA
Product category	In vitrodiagnostic device, CE
CE certified	CE
Batched	yes

### **Packaging**

Minimum order qty.	5000
Type of smallest subpackaging	bag
Piece(s) / inner box	500
Piece (s) / outer case	5000
Piece(s) / pallet	180000
Depth of case	578 mm
Width of case	368 mm
Height of case	195 mm
Case volume	0.0415 cbm
Weight of product	0.0012 kg
Weight of case	6.46 kg
EAN of inner box	4038917119757
EAN case	4038917016773

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

