



Product description

Order number	93.877.410
Product description	Storage box, hinged lid, material: PP, yellow, format: 10 x 10, for 100 collection tubes, suitable for vessels with maximum dimensions of 45 x 12 mm, 5 piece(s)/bag

Product characteristics

Closure type	hinged lid
Suitable for	vessels with maximum dimensions of 45 x 12 mm
Storage capacity	100
Format	10 x 10
Autoclavable (max. °C)	121 °C

Size

Width of product	140 mm
Height of product	55 mm
Length of product	150 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
**Storage box, hinged lid, PP, format: 10 x 10, for
100 collection tubes**

Page 2

Material & colours

Product material	Polypropylene (PP)
Colour of product	yellow
Closure material	Polypropylene (PP)
Colour of cap	yellow

Purity & certification

Product category	no medical device IVD
Batched	yes

Packaging

Minimum order qty.	20
Type of smallest subpackaging	bag
Piece(s) / inner box	5
Piece (s) / outer case	20
Piece(s) / pallet	720
Depth of case	305 mm
Width of case	279 mm
Height of case	305 mm
Case volume	0.026 cbm
Weight of product	0.1835 kg
Weight of case	4.1 kg
EAN of inner box	4038917161213
EAN case	4038917161206

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-07-18

This document was prepared by EDP support and is valid without signature.



SARSTEDT