



Product description

Order number 95.064.259
Product description Box, Styrofoam, 20-place, for micro tubes, 8 pcs./pack

Product characteristics

Storage capacity 20
Format 5 x 4

Material & colours

Product material Styrofoam
Colour of product white

Purity & certification

Product category No marking necessary
Batched Lot number

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Packaging

Minimum order qty.	16
Piece(s) / inner box	8
Piece (s) / outer case	56
Piece(s) / pallet	1792
Depth of case	401 mm
Width of case	451 mm
Height of case	280 mm
Case volume	0.0504 cbm
Weight of case	1.36 kg
EAN of inner box	4038917386883
EAN case	4038917386876

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

