



### **Product description**

Order number	95.64.981
Product description	Storage box, slip-on lid, material: cardboard, white, format: 9 x 9, for 81 collection tubes, suitable for tubes up to max. 13 mm Ø, tube height between 36-45 mm, 1 piece(s)/bag

### **Product characteristics**

Closure type	slip-on lid
Suitable for	tubes up to max. 13 mm Ø, tube height between 36-45 mm
Storage capacity	81
Format	9 x 9

### **Size**

Width of product	135 mm
Height of product	45 mm
Length of product	135 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





### **Material & colours**

Product material	cardboard
Colour of product	white
Closure material	cardboard
Colour of cap	white

### **Purity & certification**

Product category	no medical device   IVD
Batched	yes

### **Packaging**

Minimum order qty.	1
Type of smallest subpackaging	bag
Piece(s) / inner box	1
Piece (s) / outer case	120
Piece(s) / pallet	1440
Depth of case	495 mm
Width of case	419 mm
Height of case	546 mm
Case volume	0.1134 cbm
Weight of product	0.0837 kg
Weight of case	11 kg
EAN of inner box	4038917095853
EAN case	4038917007863

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

