



Product description

Order number 93.838.013
Product description Rack, material: PC, transparent, format: 12 x 4, (LxWxH): 257 x 90 x 40 mm, for 48 collection tubes, suitable for tubes, S-Monovette® 13 mm Ø, 1 piece(s)/case

Product characteristics

Suitable for tubes, S-Monovette® 13 mm Ø
Storage capacity 48
Format 12 x 4
Autoclavable (max. °C) 121 °C

Size

Width of product 90 mm
Height of product 40 mm
Length of product 257 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
Rack, PC, format: 12 x 4, suitable for tubes,
S-Monovette® 13 mm Ø

Page 2

Material & colours

Product material Polycarbonate (PC)
Colour of product transparent

Purity & certification

Product category no medical device | IVD
Batched no

Packaging

Minimum order qty. 1
Type of smallest case
subpackaging
Piece(s) / inner box 1
Piece (s) / outer case 40
Piece(s) / pallet 1200
Depth of box 92 mm
Width of box 42 mm
Height of box 259 mm
Depth of case 398 mm
Width of case 398 mm
Height of case 340 mm
Case volume 0.0539 cbm
Weight of product 0.2138 kg
Weight of case 10.29 kg
EAN of inner box 4038917121699
EAN case 4038917017947

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-07-18

This document was prepared by EDP support and is valid without signature.



SARSTEDT