



Product description

Order number	93.850
Product description	Rack, material: PC, transparent, format: 6 x 3, (LxWxH): 137 x 70 x 40 mm, for 18 collection tubes, suitable for tubes, S-Monovette® 15 mm and 13 mm Ø, 1 piece(s)/case

Product characteristics

Suitable for	tubes, S-Monovette® 15 mm and 13 mm Ø
Storage capacity	18
Format	6 x 3
Autoclavable (max. °C)	121 °C

Size

Width of product	70 mm
Height of product	40 mm
Length of product	137 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
Rack, PC, format: 6 x 3, suitable for tubes,
S-Monovette® 15 mm and 13 mm Ø

Page 2

Material & colours

Product material Polycarbonate (PC)
Colour of product transparent

Purity & certification

Product category no medical device | IVD
Batched no

Packaging

Minimum order qty. 1
Type of smallest case
subpackaging
Piece(s) / inner box 1
Piece (s) / outer case 36
Piece(s) / pallet 2376
Depth of box 72 mm
Width of box 42 mm
Height of box 138 mm
Depth of case 438 mm
Width of case 268 mm
Height of case 155 mm
Case volume 0.0182 cbm
Weight of product 0.0969 kg
Weight of case 4.55 kg
EAN of inner box 4038917084369
EAN case 4038917003223

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-07-18

This document was prepared by EDP support and is valid without signature.



SARSTEDT