



### **Product description**

Order number	93.854
Product description	Block Rack D12, for 50 collection tubes, Ø opening: 12 mm, format: 5 x 10, grey, material: PP

### **Product characteristics**

Label/ Print	with print
Number of wells	50
Storage capacity	50
Format	5 x 10
Autoclavable (max. °C)	121 °C

### **Size**

Diameter of opening	12 mm
Width of product	80 mm
Height of product	45 mm
Length of product	152 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





### **Material & colours**

Product material Polypropylene (PP)  
Colour of product grey

### **Purity & certification**

Product category no medical device | IVD  
Batched no

### **Packaging**

Minimum order qty. 1  
Type of smallest subpackaging -  
Piece(s) / inner box 1  
Piece (s) / outer case 20  
Piece(s) / pallet 1560  
Depth of case 428 mm  
Width of case 323 mm  
Height of case 130 mm  
Case volume 0.018 cbm  
Weight of product 0.1214 kg  
Weight of case 2.8 kg  
EAN of inner box 4038917084406  
EAN case 4038917003247

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

