



### **Product description**

Order number	93.870
Product description	Rack, material: PC, transparent, format: 6 x 3, (LxWxH): 137 x 70 x 40 mm, for 18 collection tubes, suitable for reaction tubes 2 ml, Microvette® , 1 piece(s)/case

### **Product characteristics**

Suitable for	reaction tubes 2 ml, Microvette®
Storage capacity	18
Format	6 x 3
Autoclavable (max. °C)	121 °C

### **Size**

Width of product	70 mm
Height of product	40 mm
Length of product	137 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Rack, PC, format: 6 x 3, suitable for reaction tubes 2 ml, Microvette®**

**Page 2**

**Material & colours**

Product material Polycarbonate (PC)  
Colour of product transparent

**Purity & certification**

Product category no medical device | IVD  
Batched no

**Packaging**

Minimum order qty. 1  
Type of smallest subpackaging case  
Piece(s) / inner box 1  
Piece (s) / outer case 36  
Piece(s) / pallet 2376  
Depth of box 72 mm  
Width of box 42 mm  
Height of box 138 mm  
Depth of case 438 mm  
Width of case 268 mm  
Height of case 155 mm  
Case volume 0.0182 cbm  
Weight of product 0.1061 kg  
Weight of case 4.88 kg  
EAN of inner box 4038917084468  
EAN case 4038917003278

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-07-18

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**