



Product description

Order number	93.893.100
Product description	Rack, material: PC, transparent, format: 10 x 2, (LxWxH): 327 x 72 x 60 mm, for 20 collection tubes, suitable for tubes 100 x 21.5 mm, 1 piece(s)/case

Product characteristics

Suitable for	tubes 100 x 21.5 mm
Storage capacity	20
Format	10 x 2
Autoclavable (max. °C)	121 °C

Size

Width of product	72 mm
Height of product	60 mm
Length of product	327 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Material & colours

Product material Polycarbonate (PC)
Colour of product transparent

Purity & certification

Product category no medical device | IVD
Batched no

Packaging

Minimum order qty. 1
Type of smallest subpackaging case
Piece(s) / inner box 1
Piece (s) / outer case 40
Piece(s) / pallet 800
Depth of case 584 mm
Width of case 387 mm
Height of case 360 mm
Case volume 0.0814 cbm
Weight of product 0.2524 kg
Weight of case 11.57 kg
EAN of inner box 4038917123136
EAN case 4038917362375

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

