



### **Product description**

Order number	92.2620.302
Product description	length: 1.5 m, material: Bioprene, transparent, for IVARO, 1 piece(s)/bag

### **Product characteristics**

Suitable for	IVARO
--------------	-------

### **Size**

Diameter	0.8 mm
Length of product	1.5 m

### **Material & colours**

Product material	Bioprene
Colour of product	transparent

### **Purity & certification**

Product category	no medical device   IVD
Batched	no

### **Packaging**

Minimum order qty.	1
Type of smallest subpackaging	bag
Piece(s) / inner box	1
Piece (s) / outer case	10
Depth of case	280 mm
Width of case	220 mm
Height of case	100 mm
Case volume	0.0062 cbm
Weight of product	0.0005 kg
EAN of inner box	4038917584197
EAN case	4038917584180

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

