



Product description

Order number	92.170.390
Product description	Transport case GS 100, for Microbial air samplers, 1 piece(s)/case

Product characteristics

Suitable for	Microbial air samplers
--------------	------------------------

Purity & certification

Product category	no medical device IVD
Batched	no

Packaging

Minimum order qty.	1
Type of smallest subpackaging	case
Piece(s) / inner box	1
Piece (s) / outer case	1
Piece(s) / pallet	1
Case volume	0.0286 cbm
Weight of product	2.7395 kg

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
Transport case GS 100, for Microbial air samplers

Page 2

Weight of case	2.74 kg
EAN of inner box	4038917101295
EAN case	4038917051590

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-07-18

This document was prepared by EDP support and is valid without signature.



SARSTEDT